

This Page Is Inserted by IFW Operations
and is not a part of the Official Record

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.


Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

IMAGES ARE BEST AVAILABLE COPY.

**As rescanning documents *will not* correct images,
please do not report the images to the
Image Problem Mailbox.**

Library of mRNA analytes



Library of mRNA analytes bound to an array





FIG. 1


(A) RNA substrate



(B) Fragmentation of RNA substrate



(C) addition of tails (UDTs) to RNA fragments



(D) Detection of tails (UDTs) on RNA fragments by binding a reagent containing signal groups (S)





FIG. 2


(A) RNA Substrate



(B) Binding of Primer to RNA Substrate



(C) Extension of Primer using RNA as template



(D) Template Independent Extension of Primer






FIG. 3


(A) RNA substrate




(B) Binding of Primer to RNA Substrate



(C) Extension of Primer using RNA as template



(D) Binding of random primers to 1st cDNA strand



(E) Extension and strand displacement of random primers






FIG. 4


(1) Binding of random primers to 1st cDNA strand



(2) Extension of random primers using 1st cDNA strand as template



(3a) Creation of functional promoters by binding of complementary strand



(3a) Creation of functional promoters by self-complementary sequences





FIG. 5


(A) Binding of Primer to analyte




(B) Extension of Primer using analyte as template



(C) Template Independent addition of dCTP



(D) Use of 3' end of 1st cNA strand for binding of Primer with T7 promoter



(E) Binding of Primer with T7 promoter to internal sequenced of cDNA





FIG. 6


(A) RNA Substrate




(B) Binding of Primer to RNA substrate



(C) Extension of Primer using RNA as template



(D) Template Independent addition of dCTP by Terminal Transferase



(E) Use of 3' end of 1st cDNA strand for binding of Primer with T7 promoter





FIG. 7


1) analyte 5' _____ A A A A A A A A 3'

2) cNA copy made from analyte 5' _____ A A A A A A A A 3'
3' _____ T T T T T T T T 5'


3a) double-stranded oligonucleotide ligated to RNA/DNA hybrid by T4 DNA ligase



3b) single-stranded oligonucleotide ligated to a single-stranded 3' tail by T4 RNA ligase




3c) double-stranded oligonucleotide ligated to single-stranded 3' tail by T4 DNA ligase




9/32


(A) RNA Substrate



(B) Binding of Primer to RNA Substrate



(C) Extension of Primer using RNA as template



(D) Nicking of cDNA strand followed by release from RNA template



(E) Template independent addition of dCTP and binding of primer with T7 Promoter






FIG. 9

10/32


(A) Capture of poly A mRNA with oligo-T beads



(B) Extension of Oligo T with poly A mRNA as template



(C) Removal of poly A mRNA and binding of random primers with T7 promoter sequence



(D) Extension of primers and strand displacement of extended primers






FIG. 10


(A) Sandwich Capture of poly A mRNA



(B) Extension of Oligo T with poly A mRNA as template



(C) Removal of poly A mRNA and binding of random primers with T7 promoter sequence



(D) Extension of primers and strand displacement of extended primers






FIG. 11


(A) Capture of poly A mRNA with oligo-T beads




(B) Extension of Oligo T with poly A mRNA as template



(C) Removal of poly A mRNA and binding of random primers to 1st cDNA strand



(D) Extension of random primers and strand displacement of extended primers



(E) Hybridization of oligo-T/T7 Pro primers to 2nd cDNA strands



(F) Extension of oligo-T/T7 Pro primers and 2nd cDNA strands





FIG. 12


(A) Poly A RNA Target




(B) Ligation of UDT to 5' end of Target




(C) Binding of primer to 3' end of Target



(D) Extension of primer



(E) Addition of Primers for Isothermal Amplification



(F) Unit length Isothermal Amplicon




FIG. 13


(A) Poly A RNA Target



(B) Ligation of UDT to 5' end of Target



(C) Binding of primer to 3' end of Target



(D) Extension of primer



(E) Addition of SDA Primer



(F) Unit length SDA Amplicon






FIG. 14


(A) Poly A RNA Target



(B) Binding of primer to 3' end of Target




(C) Extension of primer



(D) Addition of UDT (Q) to 3' end of first copy



(E) Addition of Primer for binding to Q



(F) Unit length Amplicon





FIG. 15

(1) Array with SPE's complementary to analyte "P" at Locus P, SPE's complementary to analyte "Q" at Locus Q and SPE's complementary to analyte "R" at Locus R and with UPE's comprising Poly T sequences at all three loci



(2) Binding of analyte "P" to corresponding SPE at Locus P






FIG. 16

Binding of an analyte to an array with SPE's and UPE's

(1) Analyte "P" bound to
SPE at Locus P



(2) Extension of SPE at
Locus P with analyte
"P" used as a template



(3) Removal of analyte "P"
template from extended
SPE at Locus P




FIG. 17
Extension of an SPE






FIG. 18

Binding of a UPE to an extended SPE followed by extension of the UPE

19/32

(1) Extended UPE bound to extended SPE








FIG. 19

Binding of extended SPE's and UPE's to un-extended SPE's and UPE's


*FIG. 20*

Amplification Array for Comparative Analysis

(1) Digestion of DNA with restriction enzyme R1



(2) Ligation of UDE's to DNA fragments



(3) Binding and extension of SPE primers with different 3' ends followed by extensions with UPE primers





FIG. 21

Use of an array with SPE's and UPE's for SNP analysis





FIG. 22






FIG. 23



Extension of primer P2 by
using extended SPE-P1 as
a template

Extension of primer Q2 by
using extended SPE-Q1 as
a template



Denaturation followed by annealing of
primer P2 to extended SPE-P1 and
hybridization of extended P2 to
un-extended SPE-P1

Denaturation followed by annealing of
primer Q2 to extended SPE-Q1 and
hybridization of extended Q2
to un-extended SPE-Q1

FIG. 24






FIG. 25

A) Mixture of Library of Peptide analytes with
a Library of Peptide/NAs

B) Binding of Peptide analyte #1 to Peptide/NA #1



C) Binding of Peptide/NAs to matrix through complementary sequences




FIG. 26





FIG. 27

28/32



1. Transcription product

2. 1kb PLUS DNA ladder

2000 bp
1650 bp

FIG. 28




1. Transcription product

2. 1kb PLUS DNA ladder

FIG. 29

29/32



1. Random primers -2 μ l

2. T7-C9 primers
without TdT tailing -2 μ l

3. T7-C9 primers
after TdT tailing -2 μ l

4. 1 kb PLUS DNA Ladder


5. Random primers -10 μ l

6. T7-C9 primers
without TdT tailing -10 μ l

7. T7-C9 primers
after TdT tailing -10 μ l

8. 1 kb PLUS DNA Ladder

FIG. 30




1. Taq pol. 1 cycle
2. Taq pol. 5 cycles
3. Tth pol. 1 cycle
4. Tth pol. 5 cycles
5. 1 kb PLUS DNA Ladder

FIG. 31



- | | |
|---|---|
| 1. SAMPLE 1-4 μ i TRANSCRIPTION PRODUCT | 8. SAMPLE 1-10 μ i TRANSCRIPTION PRODUCT |
| 2. SAMPLE 1-1 μ i DNA TEMPLATE | 9. SAMPLE 1-2.5 μ i DNA TEMPLATE |
| 3. SAMPLE 2-4 μ i TRANSCRIPTION PRODUCT | 10. SAMPLE 2-10 μ i TRANSCRIPTION PRODUCT |
| 4. SAMPLE 2-1 μ i DNA TEMPLATE | 11. SAMPLE 2-2.5 μ i DNA TEMPLATE |
| 5. SAMPLE 3-4 μ i TRANSCRIPTION PRODUCT | 12. SAMPLE 3-10 μ i TRANSCRIPTION PRODUCT |
| 6. SAMPLE 3-1 μ i DNA TEMPLATE | 13. SAMPLE 3-2.5 μ i DNA TEMPLATE |
| 7. 1 kb PLUS DNA LADDER | 14. 1 kb PLUS DNA LADDER |

FIG. 32



1. 1 kb PLUS DNA LADDER
2. - - - -
3. SUPERScript II (LIFE TECHNOLOGIES)
4. M-MuLV (LIFE TECHNOLOGIES)
5. M-MuMuLV (NEW ENGLAND BIOLABS)
6. ENHANCED AMV (SIGMA)
7. AMV (LIFE TECHNOLOGIES)
8. AMV (SIGMA)
9. OMNIScript (QIAGEN)
10. DISPLAY THERMO-RT (DISPLAY SYSTEMS BIOTECH)
11. POWERScript (CLONTECH)
12. - - - -
13. λ - HIND 111 MARKER